NPPTL COVID-19 Response: International Respirator Assessment

Manufacturer: 3M

Model Tested: 9502+

Date Tested: June 23, 2020

These findings pertain to the 3M, model 9502+. The packaging and labeling for this product indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Ten respirators were submitted for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found [here](#).

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 99.23% and 98.48%, respectively. All ten respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

This product has head bands/straps. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for [Crisis Capacity Strategies \(during known shortages\)](#).

Evaluation of International Respirators

Test: Modified TEB-APR-STP-0059

Date Tested: June 23, 2020

Report Prepared: June 23, 2020

Manufacturer: 3M

Item Tested: 9502+

Country of Certification: China (GB2626-2006)

Pictures have been added to the end of this report.

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
1	85	7.9	0.70	1.67	98.33
2	85	7.9	0.49	1.23	98.77
3	85	8.8	0.34	0.77	99.23
4	85	7.7	0.54	1.45	98.55
5	85	8.0	0.41	1.08	98.92
6	85	8.3	0.45	1.22	98.78
7	85	8.4	0.38	0.91	99.09
8	85	7.6	0.58	1.52	98.48
9	85	8.8	0.33	0.85	99.15
10	85	8.7	0.37	0.97	99.03
Minimum Filter Efficiency: 98.48			Maximum Filter Efficiency: 99.23		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NPPTL COVID-19 Response: International Respirator Assessment

NPPTL COVID-19 Response: International Respirator Assessment

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FITTING INSTRUCTIONS:

- Place the respirator against your face with the nosepiece across the bridge of your nose. Hold the respirator under your face and position the respirator under the chin.
- Pull the top strap over your head resting it high at the top of your head.
- Pull the bottom strap over your head position it around the neck, below the ear.
- Place the top edge of both hands at the top of the metal nosepiece. Mold the nose-piece to the shape of the nose bridge. Nosepieces using only one hand may result in a less effective respirator fit.
- Cover the front of the respirator with both hands, being careful not to disturb the position of the respirator.
 - While straps work the straps back along the side of your hand. If you CANNOT achieve a proper seal, DO NOT enter the contaminated area. See your supervisor.
 - It is very important to press the nosepiece firmly to the nose bridge to form a good seal.

WARNING

This respirator helps protect against certain particles. Misuse can result in sickness or death. For correct use, see supervisor or User Instructions on package or call 3M, 400-940-3595.

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IMPORTANT

Before use, the wearer must read and understand these User Instructions. Keep these instructions for reference. See information supplied by the manufacturer.

Use For

Particles such as those from grinding, sanding, sweeping, sawing, bagging or processing minerals, silica, coal, iron ore, heavy metals, flour, wood, pollen and certain other substances. Liquid or non-oil based particles from sprays that do not also emit oil aerosols or vapors. Metal fumes produced from welding, brazing, cutting and other operations involving heating of metal.

Do Not Use For

Do not use for gases and vapors, oil aerosols, or sandblasting, particulate concentrations that exceed either 10 times the occupational exposure limit or applicable government regulations, whichever is lower. This respirator does not supply oxygen. Do not use in situations where the oxygen concentration is less than 19.5% oxygen.

Do not use for peak-to-peak concentrations that exceed the maximum allowable concentration required by the applicable regulation and standard. The maximum allowable concentration is obtained by multiplying the Assigned Protection Factor (APF) provided by 609.138-04 Selection, Use and Maintenance of Respiratory Protective Equipment with the permissible Occupational Exposure Limit (OEL) of the contaminant.

Use Limitations and Warnings

- The wearer must first be trained in the proper use and fit test prior to using this respirator.
- Do not use this respirator when concentrations of contaminants are so high as to cause irritation to the eyes or breathing apparatus.
- Leave the contaminated area immediately if breathing becomes difficult, or dizziness or other distress occurs.
- Discard and replace respirator if it becomes damaged or breathing resistance becomes excessive.
- Inspect respirator before each use to ensure it is in good working condition. Examine all the respirator parts for signs of damage including the two straps and noseclip.
- Do not alter, wash, abuse or misuse this respirator.
- Do not use with beards or other facial hair or other conditions that prevent a good seal between the filter and the sealing edge of the respirator.
- Can be used up to the limits specified by applicable government regulations for this product.
- Failure to follow all instructions and warnings on the use of the respirator may result in injury or permanent disability.
- all times of exposure can reduce respirator effectiveness and may result in illness or permanent disability.
- Make sure the respirator is not crushed or squeezed by 21 times its force, as this may result in the deterioration and damage of the respirator. Especially pay more attention when the wearer is taking, running or conducting other activities.
- Outside forces acting upon the nose clip on the respirator may result in injury to the wearer.

Storage Conditions and Shelf Life

Shelf life of unopened product is 5 years from date of manufacture when stored at a temperature of 41°F (-20°C to +30°C) and at less than 80% relative humidity.
This respirator contains no components made from natural rubber latex.

温度范围: -20°C 至 +30°C
Temperature Range: -20°C to +30°C

最大相对湿度: <80%
Maximum Relative Humidity: <80%

保质期 (未开封产品): 自生产之日起五年
Shelf life of unopened product is 5 years from date of manufacture.

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3M China Co., Ltd.