## NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

Manufacturer: 3M

Model Tested: 1870 (TC-84A-3844)

Date Tested: July 27, 2020

Report Prepared: July 27, 2020

These findings pertain to the 3M, model 1870 (TC-84A-3844) submitted for testing; and may not be applicable to other stockpile facilities and/or under different environmental storage conditions. The maximum and minimum filter efficiency was 99.89% and 99.30%, respectively. All twenty respirators measured more than 95% efficiency. There was no designated expiration date.

NIOSH regulation sets the minimum quality and performance requirements for the approval of respirators (42 CFR 84). **NIOSH does not have requirements for shelf life or storage conditions for particulate-only APRs.** The approval holder (i.e. the entity that is granted the approval from NIOSH) is responsible for understanding how their products' design or performance may be affected by various use or storage conditions and must provide instruction for establishing the proper use, storage, and maintenance procedures for their approved products, which may include designating a shelf life. FFR or particulate filter packaging (such as the box) often includes NIOSH-approved user instructions, label information, and recommendations on shelf life. Additionally, some approval holders also disseminate recommendations related to storage and shelf life through resources such as user and web notices. The respirators tested in this study were generally not designed for long-term storage.

Based on research conducted by NIOSH and this limited testing, NIOSH does not have enough information to definitively know the level of protection that may be provided by respirators that 1) are stored for prolonged periods of times; 2) are stored under various storage conditions; or 3) have exceeded the approval holder's designated shelf life. Users of respirators that have exceeded the designated shelf life should be forewarned to avoid a false sense of confidence; these devices may not provide the same level of protection as those that have not exceeded the designated shelf life. We recommend contacting the approval holder(s) of the respirators in the stockpile with specific questions regarding the use of product beyond the manufacturer- designated shelf life.

The results provided in this letter are specific to the subset of NIOSH-approved N95s, past their designated shelf life, that were provided to NPPTL for evaluation.

These results will be added to the CDC guidance for [Stockpiled N95 Filtering Facepiece Respirators Beyond the Manufacturer-Designated Shelf Life](#).

## Evaluation of Stockpiled and Beyond Manufacturer-Designated Shelf Life N95s

**Test:** TEB-APR-STP-0059

**Manufacturer:** 3M

**Item Tested:** 1870 (TC-84A-3844)

**Expiration Date:** None Provided

**Manufacture Date:** None Provided

Filter	Flow Rate (LPM)	Initial Filter Resistance (mmH <sub>2</sub> O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency (%)
1	85	9.9	0.09	0.25	99.75
2	85	9.7	0.13	0.27	99.73
3	85	8.3	0.38	0.42	99.58
4	85	9.4	0.07	0.16	99.84
5	85	9.2	0.11	0.18	99.82
6	85	7.8	0.13	0.25	99.75
7	85	9.9	0.07	0.15	99.85
8	85	9.3	0.10	0.18	99.82
9	85	8.3	0.06	0.11	99.89
10	85	10.6	0.09	0.24	99.76
11	85	9.2	0.09	0.14	99.86
12	85	7.7	0.17	0.31	99.69
13	85	9.1	0.13	0.30	99.70
14	85	9.2	0.11	0.20	99.80
15	85	8.0	0.11	0.20	99.80
16	85	9.7	0.12	0.28	99.72
17	85	9.4	0.10	0.18	99.82
18	85	7.6	0.53	0.70	99.30
19	85	9.6	0.14	0.32	99.68
20	85	8.2	0.09	0.19	99.81
<b>Minimum Filter Efficiency: 99.30</b>			<b>Maximum Filter Efficiency: 99.89</b>		

NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

**3M** MTT-2020-318.1

**Health Care  
Particulate Respirator  
and Surgical Mask**

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LATEX

Meets CDC  
guidelines for  
*M. tuberculosis*  
exposure  
control

**FIBERGLASS FREE  
HYPOALLERGENIC  
FLUID RESISTANT**

**Appropriate for Laser Surgery**

3M  
ST. PAUL, MN, USA  
1-800-228-3957

3M  
ST. PAUL, MN, USA  
1-800-228-3957

**THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION:**

TC:	Protection <sup>1</sup>	Respirator	Cautions and Limitations <sup>2</sup>
84A-3844	N95	1870	ABCJMNOP

**1. PROTECTION**  
N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

**2. CAUTIONS AND LIMITATIONS**

- A - Not for use in atmospheres containing less than 19.5 percent oxygen.
- B - Not for use in atmospheres immediately dangerous to life or health.
- C - Do not exceed maximum use concentrations established by regulatory standards.
- J - Failure to properly use and maintain this product could result in injury or death.
- M - All approved respirators shall be selected, fitted, used, and maintained in accordance with MSHA, OSHA, and other applicable regulations.
- N - Never substitute, modify, add, or omit parts. Use only exact replacement parts in the configuration as specified by the manufacturer.
- O - Refer to users instructions, and/or maintenance manuals for information on use and maintenance of these respirators.
- P - NIOSH does not evaluate respirators for use as surgical masks.

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# NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

Fluid Resistant  
Résiste aux fluides  
Resistente a los fluidos  
Resistente a fluidos

---

FIBERGLASS FREE / LATEX FREE /  
HYPOALLERGENIC / LASER SURGERY  
EXEMPT DE FIBRE DE VERRE ET DE LATEX /  
HYPOALLERGIQUE / CHIRURGIE AU LASER  
LIBRE DE FIBRA DE VIDRIO / LIBRE DE  
LÁTEX / HIPOALERGÉNICO / CIRUGÍA LÁSER  
SEM FIBRA DE VIDRO / SEM LÁTEX /  
HIPOALERGÊNICO / CIRURGIA A LASER

---

Meets CDC guidelines for *M. tuberculosis*  
exposure control.  
Satisfait aux directives du CDC sur le  
contrôle de l'exposition à *M. tuberculosis*.  
Cumple con los lineamientos de CDC para el  
control de exposición a *tuberculosis M*.  
Atende às diretrizes CDC para controle de  
exposição ao *M. tuberculosis*.

---

Meets ASTM Fluid Resistant Challenge F1862,  
160mm Hg  
Satisfait à l'épreuve de résistance à la pénétration  
de liquides F1862, 160 mm Hg, de l'ASTM  
Cumple con los retos de resistencia de fluidos de  
ASTM F1862, 160mm Hg  
Atende às especificações de 160 mm Hg para  
Resistência a Fluido F1862 da ASTM.

**WARNING**  
This respirator helps protect against certain particulate  
contaminants but does not eliminate exposure to or the risk  
of contracting any disease or infection. Misuse may result  
in sickness or death. For proper use, see your supervisor  
or call 3M Helpline at 1-800-228-3957. In Canada,  
call 1-800-563-2921.

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3 6350

LOT 

# NPPTL COVID-19 Response: Beyond Shelf Life/Stockpiled Respirator Assessment

This respirator shape is a trademark of 3M  
Made in U.K. by 3M  
St. Paul, MN 55144-1000  
www.3M.com

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3M MEXICO S.A. DE C.V.